

# INSECTS TRAPS



## WASP ATTRACTANT®

Attractant for

*Vespula vulgaris*

*Vespula germanica*

### CHARACTERISTICS

Liquid attractant for wasps ready to use inside monitoring and capture devices for Hymenoptera. The ingredients used minimize the capture of non-target species.

It can be used indoor and outdoor.

### HOW TO USE

Use the liquid attractant within the appropriate multiple capture traps (cod. 563), both for monitoring activities and for mass trapping.

Check the traps weekly and replace the liquid every 10-15 days.

### NOTE

For correct use, follow the instructions on the label.

#### Technical data

Packaging 1L bottle in box of 6 pcs

Code 821